

NARINE BIOSUPPOSITORIES 160 mg

Pharmaceutical form: Vaginal suppositories 160mg.

Composition: Active ingredient: Lactobacillus acidophilus Er 317/402 Narine.

Description: Each suppository contains 4 billion viable lyophilized bacteria Lactobacillus acidophilus Er 317/402 Narine.

Pharmacological properties:

Narine Biosuppositories are designed to correct the vaginal microflora, effectively fight bacterial vaginosis, restore the vaginal microflora, and produce biologically active substances that are harmful to pathogens. Normally, the microflora of the vagina is represented mainly by lactobacilli, the so-called Doderlein sticks, which maintain an acidic environment in the vagina and inhibit the growth of opportunistic microorganisms.

Live lactobacilli included in the preparation have antagonistic activity against a wide range of pathogenic and opportunistic bacteria. Metabolize the glycogen of the vaginal epithelium to lactic acid, which maintains the pH of the vagina at 3.8-4.2. The lactic acid in high concentration creates unfavorable conditions for the vital activity of acid-sensitive pathogenic and opportunistic bacteria.

Indications:

Vaginal dysbiosis, including bacterial vaginosis, as an independent remedy or in complex therapy;
During and after therapy with antibiotics, cytostatics, corticosteroids, antifungal and antiviral drugs;
Subacute and chronic inflammation of the female genital organs;
After specific therapy for urinary tract infections (urethritis, cystitis...);
Prenatal preparation of pregnant women (for the prevention and treatment of vaginal dysbacteriosis);
Preparation for planned diagnostic and therapeutic gynecological operations;
While using contraception (IUD, COC, spermicides);
Hormone-dependent colpitis, senile, etc., as an independent remedy or against the background of specific hormonal therapy.

Contraindications:

Hypersensitivity, candidiasis.

Posology and method of administration:

For the purpose of prevention:

1 suppository 1 time per day for 5-7 days.

In inflammatory processes of the urogenital tract:

1 suppository 1 time per day for 5-10 days.

In case of violation of the purity of the vaginal secretion in pregnant women up to III-IV st.

1 suppository 1-2 times a day for 5-10 days or more, until the disappearance of clinical symptoms and the restoration of the purity of the vaginal secretion to I-II st.

In order to prevent purulent-septic complications:

1 suppository 1-2 times a day for 5-10 days (before the proposed operation or delivery).

Rehabilitation therapy after or during antibiotic use:

1 suppository 1 time per day for 5-10 days.

The course can be repeated for 3-4 months with an interval of 10-20 days.

Undesirable effects:

Allergic reactions.

Precautions and Warnings:

It is best to use suppositories just before your bedtime.

Pregnancy and Lactation:

Not contraindicated.

Package:

Paperboard box: 2 blisters of 5 strips.

Storage:

Store in temperature between 0°C and 20°C, and relative humidity not exceeded 65 %.

Legal status:

Without prescription.